

Watch out!

1. Many of these activities involve using a water baths at various temperatures. Make sure you are aware of what tube is placed in which water bath.
2. Be very careful when placing or removing test tubes from the boiling water bath. The tubes are very hot. Use the test tube clamps or hot gloves to remove the test tubes.
3. Some activities use several test tubes. Make sure your tubes are labeled well with their contents and with your group members' names.
4. Double check the directions to insure you are adding the correct enzyme and substrate to each tube.
5. Do not eat any of the materials in lab.
6. Report any injuries or burns to your instructor.